

POSITIONING TOOLS



STICKEMS, ADHESIVE BITEWING TABS

[Read More](#)

CODE 15720

Price: [Login to See Prices](#)

Categories: [DENTAL DISPOSABLE](#), [Disposable Products](#), [Positioning Tools](#)



ALL-BITE

[Read More](#)

CODE 409155

Price: [Login to See Prices](#)

Categories: [DENTAL DISPOSABLE](#), [Positioning Tools](#)



UNI-VERSE-ALL SENSOR POSITIONER KIT

Uni-Verse-All's unique sensor holders snap in and snap out. Use to take Anterior or Posterior/Bite Wing images. Each sensor holder adjusts "on the fly" to hold every size sensor, horizontally or vertically and are fully autoclavable.

[Read More](#)

CODE 40916

Price: [Login to See Prices](#)

Categories: [Intraoral X-ray](#), [Positioning Tools](#)

Tag: [universal sensor holder](#)



RAPID PARALLELING KIT WITH BITE WING

RAPiD is a right angle positioning device that assures correct film/cone alignment and prevents "cone cutting". It can be used for anterior, posterior, and bite-wing x-ray exams using any brand and size (3, 2, 1, and 0) of dental film.

[Read More](#)

CODE 408560

Price: [Login to See Prices](#)

Categories: [Dental Imaging Systems](#), [Positioning Tools](#)

Tags: [dental](#), [positioning kit](#), [Smart](#)



RAPID ENDODONTIC POSITIONER

Includes 1 Aiming Ring, 1 Indicator Arm, and 3 Bite Blocks.

[Read More](#)

CODE 40303

Price: [Login to See Prices](#)

Categories: [Dental Imaging Systems](#), [Positioning Tools](#)

Tags: [endodontic](#), [intraoral X-Ray](#), [positioning kit](#)



SMART POSITIONER FOR PHOSPHOR PLATE

Flow's patented SMART system makes using phosphor plates just as easy, familiar, and convenient as film.

Simply load your phosphor plate into the SMART Cassette and snap to close.

[Read More](#)

CODE 40854

Price: [Login to See Prices](#)

Categories: [CR Plates and Accessories](#), [Dental Imaging Systems](#), [Intraoral X-ray](#), [Positioning Tools](#)