

CONTRASTS

BARIUM PASTE 20ML

Barium Sulphate, radiopaque paste that provides excellent radiographic visualization of the true hoof wall. The tube contains 20mls of paste. The paste appears as an opaque (white line) on the image, when applied in a continuous bead from the coronary band to the ground surface. Due to the very thin depth of horn at the most proximal crown of the hoof wall the actual horn margin is rarely visible on even the softest penetrating views.


[Read More](#)

CODE C-GEN1010

Price: ~~\$20.00~~ \$15.00

Categories: [Contrasts](#), [Lead Letters](#), [Markers](#), [Scaling Tools](#), [Skin Markers](#)

Tags: [barium](#), [contrast](#), [equine](#), [hoof](#), [laminitis](#), [lateromedial](#)



MARKERS

BARIUM PASTE 20ML

Barium Sulphate, radiopaque paste that provides excellent radiographic visualization of the true hoof wall. The tube contains 20mls of paste. The paste appears as an opaque (white line) on the image, when applied in a continuous bead from the coronary band to the ground surface. Due to the very thin depth of horn at the most proximal crown of the hoof wall the actual horn margin is rarely visible on even the softest penetrating views.


[Read More](#)

CODE C-GEN1010

Price: ~~\$20.00~~ \$15.00

Categories: [Contrasts](#), [Lead Letters](#), [Markers](#), [Scaling Tools](#), [Skin Markers](#)

Tags: [barium](#), [contrast](#), [equine](#), [hoof](#), [laminitis](#), [lateromedial](#)



X-RAY EQUIPMENT

ATX PORTABLE X-RAY GENERATOR 5KW

ATX Portable X-ray system

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

Lead time: Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

[Read More](#)

CODE E-GEN1123

Price:

Categories: [X-Ray Components](#), [X-Ray Equipment](#)

Tags: [battery](#), [equine](#), [generator](#), [Mobile Generator](#), [Portable Generator](#), [xray](#)

ATX PORTABLE X-RAY GENERATOR 1.8KW

ATX Portable X-ray system

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

Lead time: Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

[Read More](#)

CODE E-GEN1125

Price:

Categories: [X-Ray Components](#), [X-Ray Equipment](#)

Tags: [battery](#), [equine](#), [generator](#), [Mobile Generator](#), [Portable Generator](#), [xray](#)

ATX PORTABLE X-RAY GENERATOR 2.4KW

ATX Portable X-ray system

- Designed with innovation in mind, Ultra-light, high frequency x-ray system.
- Compact size and lightweight body
- Durable rugged exterior

Provided with Compliance testing as per local radiation health regulations.

Lead time: Subject to stock availability, transit times and radiation licensing. Estimated shipping date to be provided once the purchase is confirmed

[Please contact us for pricing](#)

[Read More](#)

CODE E-GEN1124

Price:

Categories: [X-Ray Components](#), [X-Ray Equipment](#)

Tags: [battery](#), [equine](#), [generator](#), [Mobile Generator](#), [Portable Generator](#), [xray](#)



