

MARKERS

BEEKLEY CT-SPOT® 119

BEEKLEY CT-SPOT® 119

Ref: 119, (QTY: 85/box)

Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.

- 2.3mm non-metallic pellet
- for 3-points set-ups and isocenters

[Read More](#)

CODE C-BEEK119

Price: \$183.00

Categories: [Markers](#), [Skin Markers](#)

Tags: [CT Simulation Skin Markers](#), [Grids and Markers](#)



