

# DENTAL

## UNI-VERSE-ALL SENSOR POSITIONER KIT

Uni-Verse-All's unique sensor holders snap in and snap out. Use to take Anterior or Posterior/Bite Wing images. Each sensor holder adjusts "on the fly" to hold every size sensor, horizontally or vertically and are fully autoclavable.

[Read More](#)

**CODE** C-DEN1027

**Price:** \$95.00

**Categories:** [Dental](#), [Dental Positioning Aids](#)

**Tag:** [universal sensor holder](#)



