

UNCATEGORISED



VETARIO INTENSIVE CARE T50M

- Overall height (mm) 560
- Overall width (mm) 690
- Overall depth (mm) 490
- Floor area (mm) 600 x 400
- Effective volume (litres) 100
- Power consumption (Watts) Max 160, Typical 95
- Power supply 230v, 50Hz
- Weight (Kg) 9.2

[Read More](#)
CODE 2VetT50M
Price: ~~\$2,650.00~~ \$2,600.00
Category: [Uncategorised](#)
Tags: [incubator](#), [T50M](#), [vetario](#), [vetariot50m](#)

