

POSITIONING AIDS

RADIOGRAPHY CALIPER

Easy to use caliper for anatomy thickness measurements.

Set you cm thickness generator factors for optimum plate dose.


[Read More](#)

CODE E-ACC1033

Price: \$125.00

Categories: [Foams](#), [Positioning Aids](#)

Tags: [anatomy](#), [factors](#), [measurement](#), [radiography](#), [thickness](#), [X-ray](#)



X-RAY EQUIPMENT

DENTAL X-RAY PACKAGE CR2

The Veterinary Dental X-ray Package is comprised of:

- Owandy RX Intra-oral Wall mount x-ray [View here](#)
- APIXIA EXM [View here](#)
- or APIXIA EXL [View here](#)
- Initial 1 year warranty & unlimited remote support

Added benefits:

- Lead gown start up package (1 lead gown, 1 thyroid collar, 1 pair of vet mitts)

[Read More](#)

CODE N/A

Price:

Categories: [X-Ray Equipment](#), [X-Ray Full Solutions](#)

Tags: [Cr](#), [CR2](#), [dental](#), [X-ray](#), [X-ray room](#), [Xray room](#)

OWANDY OPTEO DR SENSORS (USB)

Designed to stay with you for the long term, the Opteo digital x-ray sensor is sturdy and resistant. The cable connecting it to the USB interface unit has been specially reinforced to offer greater resistance to impacts and twisting.

Easy to use, Opteo offers the best image quality/price ratio. High-quality images with minimal "noise" make for a fast, exact diagnosis every time.

Size 1: 25,2mm x 36,9mm

Size 2: 30,8mm x 43,2mm

Sensors can be supplied with a computer at extra cost. Please contact ATX on **1800 330 118** if you would this option.

[Read More](#)

CODE N/A

Price:

Categories: [Digital Radiography](#), [X-Ray Equipment](#)

Tags: [Dental Sensor](#), [Dental xray](#), [Opteo](#), [Opteo DR](#), [X-ray](#)

DENTAL X-RAY DR SOLUTION

The ATX DENTAL X-RAY DR PACKAGE is the perfect choice for in-practice dental x-raying, Owandy RX has been constructed for Accuracy and precision. Highly precise exposure time, accurately controlled through a timer microchip. Safeguard your patients with HF technology, which helps to reduce x-ray doses by eliminating harmful soft rays. Benefits: Precision, Reliability, Ergonomic design, intuitive and easy to install

Lead time: Subject to stock availability, transit times and radiation licensing . Estimated shipping date to be provided once the purchase is confirmed.

[Read More](#)

CODE N/A

Price:

Categories: [X-Ray Equipment](#), [X-Ray Full Solutions](#)

Tags: [dental](#), [X-ray](#), [X-ray room](#), [Xray room](#)



VETERINARY MOBILE DENTAL X-RAY PACKAGE

The Veterinary Dental Mobile X-ray Package is comprised of:

- REMEX handheld X-ray generator [View here](#)
- Owandy Opteo HD DR Digital detector (USB) [View here](#)
- Initial 1 year warranty & unlimited remote support
- ATX Cloud PACS [View here](#)

Added benefits:

- Lead gown start up package (1 lead gown, 1 thyroid collar, 1 pair of vet mitts)

[Please contact us for pricing](#)

[Read More](#)

CODE N/A

Price:

Categories: [X-Ray Equipment](#), [X-Ray Full Solutions](#)

Tags: [dental](#), [mobile xray](#), [X-ray](#), [X-ray room](#), [Xray room](#)



