

# UNCATEGORISED



## DENTAL X-RAY TRAINING PHANTOM

## DENTAL X-RAY TRAINING PHANTOM

Dental X-Ray Training Procedure Phantom

The dental intra-oral X-ray training phantom includes a special DENTOFORM® model.

The dental intra-oral X-ray training phantom includes a special DENTOFORM® model with radio-opaque metal teeth, flexible finger for holding film, bite-opening instruments, and latex tongue.

Adult version. Paediatric version also available.

### Technical Specifications

- FEATURES: Chair mount included
- WEIGHT: 3 kg
- OPTIONAL: Holder for bench/table

[Read More](#)

**CODE** 12704

**Price:** \$10,000.00

**Category:** [Uncategorised](#)

