

UNCATEGORISED



DENTAL X-RAY TRAINING PHANTOM

DENTAL X-RAY TRAINING PHANTOM

Dental X-Ray Training Procedure Phantom

The dental intra-oral X-ray training phantom includes a special DENTOFORM® model.

The dental intra-oral X-ray training phantom includes a special DENTOFORM® model with radio-opaque metal teeth, flexible finger for holding film, bite-opening instruments, and latex tongue.

Adult version. Paediatric version also available.

Technical Specifications

- FEATURES: Chair mount included
- WEIGHT: 3 kg
- OPTIONAL: Holder for bench/table

[Read More](#)

CODE 12704

Price: \$10,000.00

Category: [Uncategorised](#)

