

# UNCATEGORIZED

## OPG BITE BLOCKS

OPG Bite Block to fit. Sold individually

[Read More](#)

**CODE** C-DEN1014

**Price:** \$1.65

**Category:** [Uncategorized](#)

**Tag:** [OPG & Dental Disposable](#)



# CBCT 3D

3D imaging for dental, periodontic, orthodontic and maxillofacial.

## CBCT 3D IMAGING

Rotograph Evo 3D is a new dental imaging system with 3-in-1 functionality: Panoramic, Cephalometric, 3D.

This system is based on the highly regarded mechanical platform of the Rotograph Evo, adapted mechanically and electronically to operate with the latest "Cone Beam" technology allowing 3D reconstruction of the jaws.

Core of the Rotograph Evo 3D is the innovative Flat Panel Detector (FPD) using Amorphous Silicon technology associated with Cesium Iodide (CsI) scintillator.

Thanks to its peculiarity to work in either Panoramic or 3D modality, Pan images are acquired directly and not reconstructed from a volume, while the unit's 3D Field of View of 8,5×8,5 cm allows to obtain the entire dental status of the patient in a single exposure. The

Rotograph Evo 3D can also be integrated with a digital Cephalometric arm, or upgraded at a later stage in case of need, in order to expand applications to the orthodontic field.

[Read More](#)

**CODE** N/A

**Price:** \$125,000.00

**Categories:** [CBCT 3D](#), [Dental Imaging Systems](#)



# DENTAL

## VIDEOGRAPH

Leveraging an enlarged experience in all sectors of dental radiology, including panoramic and 3D cone beam technologies, Villa S.M. includes in its product portfolio the Videograph, a hi-tech intraoral sensor available in two different models: regular and HD.

[Read More](#)

**CODE** 9458800021

**Price:** \$8,500.00

**Category:** [Dental Imaging Systems](#)



## OPG EAR OLIVE

OPG Ear Olive, Pack of 50.

[Read More](#)

**CODE** C-DEN1015

**Price:** \$130.00

**Category:** [Dental](#)

**Tag:** [OPG Disposable](#)



# DENTAL EQUIPMENT



## DIGITAL ACQUISITION - CR

[Read More](#)

**CODE** N/A

**Price:** \$9,500.00

**Categories:** [Dental Equipment](#), [Dental Imaging Systems](#), [Intraoral X-ray](#)



# DENTAL IMAGING SYSTEMS

## DIGITAL ACQUISITION - CR

[Read More](#)

**CODE** N/A

**Price:** \$9,500.00

**Categories:** [Dental Equipment](#), [Dental Imaging Systems](#), [Intraoral X-ray](#)



## CBCT 3D IMAGING

Rotograph Evo 3D is a new dental imaging system with 3-in-1 functionality: Panoramic, Cephalometric, 3D.

This system is based on the highly regarded mechanical platform of the Rotograph Evo, adapted mechanically and electronically to operate with the latest "Cone Beam" technology allowing 3D reconstruction of the jaws.

Core of the Rotograph Evo 3D is the innovative Flat Panel Detector (FPD) using Amorphous Silicon technology associated with Cesium Iodide (CsI) scintillator.

Thanks to its peculiarity to work in either Panoramic or 3D modality, Panoramic images are acquired directly and not reconstructed from a volume, while the unit's 3D Field of View of 8,5×8,5 cm allows to obtain the entire dental status of the patient in a single exposure. The

Rotograph Evo 3D can also be integrated with a digital Cephalometric arm, or upgraded at a later stage in case of need, in order to expand applications to the orthodontic field.

[Read More](#)


**CODE** N/A

**Price:** \$125,000.00



**Categories:** [CBCT 3D](#), [Dental Imaging Systems](#)

# INFECTION CONTROL



## OPG CENTERING POST

75x6mm OPG Centering Post. Sold individually

[Read More](#)

**CODE** C-DEN1013

**Price:** \$2.00

**Categories:** [Disposable Products](#), [Infection Control](#), [OPG & Dental Disposable](#)

**Tag:** [OPG Centering Post](#)

# INTRAORAL X-RAY

## DIGITAL ACQUISITION - DR

[Read More](#)

**CODE** N/A

**Price:** \$9,500.00

**Category:** [Intraoral X-ray](#)



## DIGITAL ACQUISITION - CR

[Read More](#)

**CODE** N/A


**Price:** \$9,500.00

**Categories:** [Dental Equipment](#), [Dental Imaging Systems](#), [Intraoral X-ray](#)





# OPG & DENTAL DISPOSABLE



## OPG CENTERING POST

75x6mm OPG Centering Post. Sold individually

[Read More](#)

**CODE** C-DEN1013

**Price:** \$2.00

**Categories:** [Disposable Products](#), [Infection Control](#), [OPG & Dental Disposable](#)

**Tag:** [OPG Centering Post](#)

# X-RAY EQUIPMENT

## ENDOGRAPH DC MOBILE VERSION

Endograph DC matches any practice, as it is available with 3 different sizes of extension arms that can be adapted to virtually any room. Upon request, Endograph DC can come also as a mobile unit mounted on trolley, particularly indicated in emergency contexts. Additionally, the standard availability of preset settings for films, Photostimulable Phosphor Plates (PSP), or Digital Detectors, allows the use with any receptor, from any manufacturer.

[Read More](#)

**CODE** 9459000313

**Price:** \$7,000.00

**Category:** [X-Ray Components](#)

**Tags:** [intraoral X-Ray](#), [Mobile](#), [wireless](#)



## OWANDY RX INTRAORAL X-RAY

The Owandy RX has been constructed for Accuracy and precision. Highly precise exposure time, accurately controlled through a timer microchip. Safeguard your patients with HF technology, which helps to reduce x-ray doses by eliminating harmful soft rays.

[View OWANDY RX Brochure](#)

[Please contact us for pricing](#)

[Read More](#)

**CODE** N/A

**Price:**

**Categories:** [X-Ray Components](#), [X-Ray Equipment](#)

**Tags:** [intraoral](#), [Owandy](#), [RX Intraoral](#), [xray](#)



